



國立清華大學  
NATIONAL TSING HUA UNIVERSITY


**daw**i Data Analysis &  
Interpretation Lab

# Readme – motion correction

*registration\_ver4*


**Edit by Wun-Ci Chen, Kai-Chun Jhan**

**2024/06/21**




# 1. System requirements

- Our environment
  - Windows 10
  - Matlab R2021b
  - CPU 11<sup>th</sup> Gen Intel Core
- This code has been tested on our environment and Matlab R2020b.
- No non-standard software required.




## 2. Installation guide

- Please refer to the following URL for MATLAB installation instructions.
  - <https://www.mathworks.com/help/install/install-products.html>
- Typically, installing MATLAB on a "normal" desktop computer takes around 30 minutes to 1 hour.
- Typically, installing our code takes less than one minute.




# 3. Demo

- Perform motion correction using auxiliary corrected data (denoised, cropped, etc.), then apply the parameters back to the original data.
- Running the program with demo data on a normal computer should complete within five minutes.




# 3. Demo




We use the DeepCAD-RT algorithm to obtain denoised images.  
<https://github.com/cabooster/DeepCAD-RT>

## 4. Instructions for use











# Folder organization



# Example



Run (user defined)

 img	2024/6/21 下午 01:53	檔案資料夾	
 img_trial	2024/6/21 下午 01:43	檔案資料夾	
 aligntrials	2024/6/21 下午 02:16	MATLAB Code	3 KB
 apply_raw	2024/6/21 下午 02:09	MATLAB Code	2 KB
 registration_trial	2024/6/21 下午 02:09	MATLAB Code	3 KB
 run_loop_register	2024/6/21 下午 02:09	MATLAB Code	2 KB
 run_loop_register_all	2023/3/6 下午 01:59	MATLAB Code	1 KB
 transfer_stack_tif	2023/3/6 下午 02:30	MATLAB Code	1 KB

Store all the code and create two folders named "img" and "img\_trial".



# Example




2024/6/21 下午 02:15

檔案資料夾


2024/6/21 下午 02:15

檔案資料夾

Under the newly created "img" folder, create folders named after each data ID  
(note: naming order may vary).




# Example



ABCDE (user defined)


- Inside each group's folder, create 2 folders as shown in the attached image.
- In the "img\_raw" folder, place raw data for this group (multiple channels can be placed at once).
- In the "img\_auxiliary" folder, place auxiliary corrected data for this group.
- Data format must be in 3D stack (x-y-t).

\* Multiple channels must be acquired simultaneously with identical displacements in time.




- Go back to the top-level "img\_trial" folder created earlier, and place the first registered TIFF image (or any cropped image that aids in mutual alignment of each group) for each registered group inside.
- Note: The TIFF files placed must be named the same as the folders in the third level (i.e., the names of each group's folder).

# Program description



# Folder arrangement for all

- Folder structure "run" (users can change the name):
  - Place all program files ("p files") here.
- run/img
  - Place all subfolders containing images to be processed under "img".
- run/img/ABCDE/img\_raw(or img\_auxiliary)
  - Organize images for processing under either "img\_raw" or "img\_auxiliary".
- run/img/ABCDE/singletrial\_raw(or singletrial\_auxiliary)
  - Results from processing a single group will correspondingly be stored here.
- run/img/ABCDE/aligntrial\_raw(or aligntrial\_auxiliary)
  - Final aligned results from all groups will be stored here.
- run/img\_trial
  - Copy the first TIFF file (or cropped image) from each registered group here, for alignment purposes.



# Execution method and workflow

## 1. `run_loop_register_all(false)`:

- Execute this function to perform registration for each single group.
- This will generate registration results in each group's "singletrial" folder.

The parameter for executing ``run_loop_register_all``:

- True: Adjust image contrast before applying registration parameters.
- False: Apply motion correction without adjusting image contrast.

## 2. `aligntrials`:

- Run this step to align all groups.
- It aligns all groups based on the first trial's registration results.

## 3. `registration_trial`:

- Execute this function in MATLAB. MATLAB will prompt for which trial to adjust ("Which trial?").
- Input a number equal to or greater than 2 to specify which trial to adjust.
- The first trial in alphabetical order serves as the baseline for alignment (1).
- Subsequent trials are numbered sequentially starting from 2.